

Test Case Planning

Group members	Ryan Brady, Aaron Matthews, Victoria Castagnola, Nicholas Augustus, Jacob Alexander, Elijah Borum
----------------------	---

1. Use Case 1 - *User Search for Patients*

Test Case ID	TC01 - Full Name Search
Test Objective	Patient search by exact patient name
Preconditions	Patient database has been initialized with patient data
Test steps	<ol style="list-style-type: none">1. Select search option on home screen2. Input patient's first name3. Input patient's last name4. Press the button to initiate search
Input Values	Patient's first and last name
Expected results	Patient file

Test Case ID	TC02 - Partial Name Search
Test Objective	Patient search by partial patient name
Preconditions	Patient database has been initialized with patient data
Test steps	<ol style="list-style-type: none">1. Select search option on home screen2. Input either the patient's first OR last name3. Press the button to initiate search
Input Values	Patient's first name OR patient's last name
Expected results	List of patients with either the same first or last name as the input

Test Case ID	TC03 - Search by ID
Test Objective	Patient search by Medical Record Number
Preconditions	Patient database has been initialized with patient data
Test steps	<ol style="list-style-type: none">1. Select search option on home screen2. Select the advanced search drop down and select Medical Record Number3. Enter the patient's Medical Record Number4. Press the button to initiate search
Input Values	Patient's Medical Record Number

Expected results	The file of the patient with the given number
-------------------------	---

Test Case ID	TC04 – Search by Room Number
Test Objective	Patient search by Room Number
Preconditions	Patient database has been initialized with patient data
Test steps	<ol style="list-style-type: none">1. Select search option on home screen2. Select the advanced search drop down and select Room number3. Enter the patient's Room Number4. Press the button to initialte search
Input Values	Patient's Room Number
Expected results	The file of the patient with the given room number